

April 2008

New Laccase Test Kit By Dolmar - Works In the Presence Of SO₂!!!

Fleurieu Winery Supplies distributes in Australia & NZ, which is an improvement on anything currently available. See the attached data sheet. Its easily posted over to you in NZ.

The Dolmar kit still uses syringaldazine but is a special formulation with other components to confer it the following advantages being :

1. The sensitivity of the test: the test will detect potential laccase oxydative activity that the other test kits do not detect. And so, the winemaker can act much more quickly/earlier with higher SO₂, addition of tannins for their antiseptic action, faster processing/winemaking decisions and thus avoid oxidation problems. 10 to 100 x more sensitive
2. The 2nd major advantage : the new kit works in presence of SO₂ (the other test kits do not) so we can test must samples even if the grapes have been treated in the bins with PMS and/or ascorbic.

For Further Information contact-

Olga Kostic

Fleurieu Winery Supplies

PO Box 7, Yankalilla 5203

Tel/Fax: +61-8-8558 3656

Mobile: +61 419 039 167

okostic@iprimus.com.au

www.newworldwinemaker.com (Anchor Wine Yeast)

www.fordras.com/www.lysoclub.com (Lysozyme)

www.ADD-FOOD.com (Natuzym enzymes)

www.oenoprod.com (Tannins)

SNIFF'Brett Brettanomyces testing kits

PROSTAB protein stability testing in 10mins

SOFA Fermentation simulation software, INRA France

ENODEFECTO the sensory taints/faults kit